

Halton and St Helens Wheelchair Service

How to ensure your battery lives a long and healthy life...

1. Charge your batteries for at least 18 hours prior to their first use.
2. Make sure you charge your battery every day, even after a small amount of use.
3. Charge your battery for as long and as often as you possibly can
4. Don't worry! Charging your batteries overnight is fine - you cannot overcharge batteries.
5. Do not charge your batteries in conditions under 10°C or over 30°C.
6. Do not store your powered wheelchair in outdoor conditions where adverse weather can affect the electrical components and batteries.
7. If possible, charge your battery in a different room to where you are sleeping (in a ventilated room)
8. Try to allow a complete charging cycle of 12 hours at least once a week.
9. Always try to charge your battery before the red light appears on your joystick controller.
10. Your battery will reach full capacity after 10 to 20 charging cycles.

And finally...

Please read your User Manual very carefully to find out other ways to keep your chair in top condition.

Halton & St Helens Wheelchair Service: 01928 582939

Ross Care Repair Service: 0151 653 6000